

Revised diagnostic reference levels in CT imaging: Progress in radiation protection?

Dr. Christoph Aberle, Medical Physicist

KSR Webinar
Radiation Protection in Medicine

Friday, 29th January 2021

2 projects were funded by the Federal Office of Public Health (FOPH).

People directly involved:

R. Treier¹, S. Schindera², F. Verdun³, A. Viry³, D. Racine³,
N. Ryckx³, T. Lima^{3,4}, C. Aberle⁵

1 Radiological Protection Division, Federal Office of Public Health, Bern, Switzerland


2 Department of Radiology, Kantonsspital Aarau, Aarau, Switzerland

3 Institute of Radiation Physics, Lausanne University Hospital and University of Lausanne, Lausanne, Switzerland

4 Department of Radiology and Nuclear Medicine, Kantonsspital Luzern, Luzern, Switzerland

5 Department of Radiology, University Hospital Basel, University of Basel, Basel, Switzerland






- Introduction
 - Radiation exposure due to CT
 - Diagnostic Reference Levels (DRLs)
- Current Swiss DRLs for CT
- Comparison to previous Swiss DRLs
- Dose and image quality
- Summary




Introduction

- Sharp increase in the number of CT scanners from 1998 to 2017
- Corresponding increase in the number of CT examinations

CT scanners in Switzerland




CT examinations per 1'000 inhabitants



Introduction

- As of 2018, CT contributes **69.4%** to the radiation doses used in radiodiagnostic procedures.
- 0.96 mSv per inhabitant per year due to CT



Introduction

- DRLs are a tool for dose optimization
- The DRL concept is implemented in the RPO

Art. 35 Diagnostic reference levels

¹ The FOPH shall publish recommendations on radiation doses for diagnostic, interventional or nuclear medicine examinations in the form of diagnostic reference levels.

² To this end, it shall conduct national surveys based on the data specified in Article 34 paragraph 2, take international recommendations into account and publish the results.

³ Licence holders must regularly review their own practices and account for any deviations from diagnostic reference levels.




Introduction

Setting up DRLs:


1. Dose distributions (CTDI_{vol} , DLP) for different examination types are obtained
2. DRLs are set at the 75th percentile

Application of DRLs:

1. If a DRL is consistently exceeded, investigations have to be performed
2. Dose levels have to be reduced if there is no good reason for elevated doses



Current Swiss DRLs for CT




- published by the FOPH in 2018 [7,8]
- 14 participating institutes
- 50 CT scanners
- 5 dose management software tools
- Data from 2014-2017 included

Institute	Number of CTs	Dose Management Software
Centre Hospitalier Universitaire Vaudois (CHUV), Lausanne	5	DoseWatch (GE)
Kantonsspital Frauenfeld (KSF), Frauenfeld	2	Radimetrics (Bayer)
Kantonsspital Graubünden (KSGR), Chur	3	Radimetrics (Bayer)
Universitätsspital Zürich (USZ), Zürich	3	tqm DOSE (Agfa)
Kantonsspital Münsterlingen (KSM), Münsterlingen	2	Radimetrics (Bayer)
Universitätsspital Basel (USB), Basel	3	Radimetrics (Bayer)
Hôpitaux Universitaires de Genève (HUG), Genève	7	Radimetrics (Bayer)
Affidea Schweiz, 5 sites	5	DoseWatch (GE)
Privatklinikgruppe Hirslanden, 9 sites	10	DoseM (INFINITT)
Stadtspital Triemli, Zürich	2	DoseWatch (GE)
Kantonsspital Aarau (KSA), Aarau	3	in-house solution
Merian Iselin Klinik, Basel	1	Radimetrics (Bayer)
Kantonsspital Baselland (KSBL), 3 sites	3	Radimetrics (Bayer)
Felix-Platter-Spital (FPS), Basel	1	Radimetrics (Bayer)

data from [8]


Current Swiss DRLs for CT



Anatomical region	Master protocol (depending on clinical indication)	
Head	Head	Head phantom (16 cm diameter)
	Facial Bones / Sinuses	
	Low Dose Sinuses (Sinusitis)	
Neck	Neck	Body phantom (32 cm diameter)
	CT Angiography (Carotid Angio)	
Chest	Chest	
	CT Angiography (Exclusion Pulmonary Embolus -PE)	
Abdomen-Pelvis	Abdomen-Pelvis	
	Exclusion Kidney Stones	
	Liver Multiphase (e.g. HCC)	
	CT Angiography (Abdomen-Pelvis)	
Chest-Abdomen-Pelvis	Chest-Abdomen-Pelvis	
	CT Angiography (Chest-Abdomen-Pelvis)	
Musculoskeletal	Cervical Spine	
	Thoracic Spine / Lumbar Spine	

data from [8]

Current Swiss DRLs for CT




Master protocol	Number of exams	Number of CTs with at least 20 exams
Head	58,323	45
Facial Bones / Sinuses	6,322	31
Low Dose Sinuses (Sinusitis)	2,328	15
Neck	3,079	30
CT Angiography (Carotid Angio)	4,528	19
Chest	32,437	49
CT Angiography (Exclusion PE)	16,998	41
Abdomen-Pelvis	33,895	49
Exclusion Kidney Stones	12,625	35
Liver Multiphase (e.g. HCC)	3,621	26
CT Angiography (Abdomen-Pelvis)	1,102	18
Chest-Abdomen-Pelvis	25,613	35
CT Angiography (Chest-Abdomen-Pelvis)	2,480	22
Cervical Spine	8,710	36
Thoracic Spine / Lumbar Spine	8,208	35
Total	220,269	data from [8]


Current Swiss DRLs for CT

- Example: abdomen-pelvis master protocol
- DRL calculated according to ICRP publication 135 [9]: 75th percentile of the distribution of the CT scanner medians

DRL: 75th percentile
50th percentile
25th percentile



Current Swiss DRLs for CT




Master protocol	25th Percentile		50th Percentile		75th Percentile (DRL)	
	CTDI _{vol} [mGy]	DLP [mGy cm]	CTDI _{vol} [mGy]	DLP [mGy cm]	CTDI _{vol} [mGy]	DLP [mGy cm]
Head	35	650	42	750	51	890
Facial Bones / Sinuses	9	150	15	240	25	420
Low Dose Sinuses (Sinusitis)	3	40	5	70	6	90
Neck	9	240	12	290	16	410
CT Angiography (Carotid Angio)	5	190	8	260	11	360
Chest	5	170	6	210	7	250
CT Angiography (Exclusion PE)	4	150	6	200	8	300
Abdomen-Pelvis	8	390	10	470	11	540
Exclusion Kidney Stones	4	150	4	180	6	280
Liver Multiphase (e.g. HCC)	7	220 [830]	9	300 [960]	11	350 [1170]
CT Angiography (Abdomen-Pelvis)	7	340	9	450	11	530
Chest-Abdomen-Pelvis	7	500	9	610	11	740
CT Angiography (Chest-Abdomen-Pelvis)	5	340	6	450	10	730
Cervical Spine	12	260	14	300	17	360
Thoracic Spine / Lumbar Spine	13	–	18	–	25	–

data from [8]

Values are per acquisition, only for liver multiphase a value for total examination
DLP is quoted in addition [in square brackets]

DRLs

Comparison to previous Swiss DRLs



The average DRL reduction is -30% for CTDI_{vol} (0% to -47%) and -22% for DLP (+20% to -40%).


Master protocol	current Swiss DRLs [8]		previous Swiss DRLs (2010) [10]	
	CTDI _{vol} [mGy]	DLP [mGy cm]	CTDI _{vol} [mGy]	DLP [mGy cm]
Head	51 (78%)	890 (89%)	65	1000
Facial Bones / Sinuses	25 (100%)	420 (120%)	25	350
Low Dose Sinuses (Sinusitis)	6	90	–	–
Neck	16 (53%)	410 (68%)	30	600
CT Angiography (Carotid Angio)	11 (55%)	360 (72%)	20	500
Chest	7 (70%)	250 (63%)	10	400
CT Angiography (Exclusion PE)	8 (53%)	300 (67%)	15	450
Abdomen-Pelvis	11 (73%)	540 (83%)	15	650
Exclusion Kidney Stones	6	280	–	–
Liver Multiphase (e.g. HCC)	11	350 [1170]	–	–
CT Angiography (Abdomen-Pelvis)	11 (73%)	530 (82%)	15	650
Chest-Abdomen-Pelvis	11 (73%)	740 (74%)	15	1000
CT Angiography (Chest-Abdomen-Pelvis)	10	730	–	–
Cervical Spine	17 (57%)	360 (60%)	30	600
Thoracic Spine / Lumbar Spine	25 (83%)	-	30	–

Round brackets: values relative to the previous DRLs

data from [8,10]

Dose and image quality

- DRLs don't measure image quality
- There is the possibility of excessive dose reduction
- If local dose values are below the 25th percentile, image quality should be considered the first priority for further optimization
- Optimization might also mean increasing the dose in some cases
- The goal is standardization of image quality




Dose and image quality

Follow-up project to study the link of DRLs to image quality:

- mathematical model observer phantom study
- The graph shows an example for low contrast lesion detection (5 mm diameter spheres with 20 HU contrast in a phantom)

→ Image quality is not very high around the DRL and it does not reach a plateau in the clinically relevant dose range

→ There is a need for specific dose optimization for each CT scanner



[11,12]

Summary

- Swiss DRLs for adult CT were updated in 2018
- There was a substantial average reduction of DRLs compared to 2010 due to new CT technology and optimization efforts.
- Average DRL reduction compared to 2010 by 30% for CTDIvol and 22% for DLP
- Dose optimization is not always dose reduction, image quality has to be standardized



References

1. Federal Office of Public Health, Switzerland, Radiation Protection Division (2017) Strahlenschutz und Überwachung der Radioaktivität in der Schweiz: Ergebnisse 2017. available via <https://www.bag.admin.ch/bag/de/home/das-bag/publikationen/taetigkeitsberichte/jahresberichte-strahlenschutz-umweltradioaktivitaet-und-dosimetrie.html>. Accessed 22 Jan 2021
2. Aroua A, Burnand B, Decka I, Vader JP, Valley JF (2002) Nationwide survey on radiation doses in diagnostic and interventional radiology in Switzerland in 1998. Health Phys 83:46–55
3. Aroua A, Trueb P, Vader JP, Valley JF, Verdun FR (2007) Exposure of the Swiss population by radiodiagnostics: 2003 review. Health Phys 92:442–448
4. Samara ET, Aroua A, Bochud FO et al (2012) Exposure of the Swiss population by medical X-rays: 2008 review. Health Phys 102:263–270
5. Le Coultre R, Bize J, Champendal M et al (2016) Exposure of the Swiss population by radiodiagnostics: 2013 review. Radiat Prot Dosimetry 169:221–224
6. Bize J, LeCoultre R, Viry A, Verdun FR (2020) Strahlenexposition der Schweizer Bevölkerung durch Röntgendiagnostik 2018. available via <https://www.bag.admin.ch/bag/de/home/gesund-leben/umwelt-und-gesundheit/strahlung-radioaktivitaet-schall/strahlenanwendungen-in-der-medizin/diagnostische-strahlenexposition-in-der-medizin.html>. Accessed 22 Jan 2021
7. Federal Office of Public Health, Switzerland, Radiation Protection Division (2018) Directive R-06–06: Diagnostische Referenzwerte in der Computertomografie. Available via <https://www.bag.admin.ch/dam/bag/de/dokumente/str/std/bag-wegleitungen/2018-06-15-r-06-06.pdf.download.pdf/R-06-06wd.pdf>. Accessed 22 Jan 2021
8. Aberle C, Ryckx N, Treier T, Schindera ST (2020) Update of national diagnostic reference levels for adult CT in Switzerland and assessment of radiation dose reduction since 2010. Eur Radiol 30:1690–1700
9. Vañó E, Miller DL, Martin CJ et al (2017) ICRP Publication 135: Diagnostic reference levels in medical imaging. Ann ICRP 46(1):1–144
10. Treier R, Aroua A, Verdun FR, Samara E, Stuessi A, Trueb PR (2010) Patient doses in CT examinations in Switzerland: implementation of national diagnostic reference levels. Radiat Prot Dosimetry 142:244–254
11. Viry A, Aberle C, Lima T, Treier R, Schindera ST, Verdun FR, Racine D (2021) Assessment of task-based image quality for abdominal CT protocols linked with national diagnostic reference levels. manuscript under review.
12. Racine D, Viry A, Aberle C, Lima T, Treier R, Schindera ST and Verdun FR: Evaluation of the image quality in abdominal CT protocols based on Swiss diagnostic reference levels. Optimisation in X-ray and Molecular Imaging 2020, O2-5.